

# EnviroMed-Bag® DELUXE

**LDI**medical™  
Sustainable Technology · Clinical Innovation · Ultimate Performance

## Eco-Friendly Body Bags

PART #	SIZE	PACKAGING
FBB2-3692H	36" x 92" Infant, curved zipper with tags	6 bags per case
FBB2-3692	36" x 92" Pediatric, curved zipper with tags	12 bags per case



### Product Details—Body Bags With or Without Handles

- **Scrim Reinforced 14 Mil Thickness Construction**—Provides superior puncture and tear resistance
- **Heavy Duty, Rust Proof Zipper**—With two pulls for easy open and close
- **Handles**—Available with or without six integral handles
- **Envelope Style Opening**—For easier insertion and removal of cadavers

### Manage Risk and Comply with OSHA Standards

- **Compliance**—Aides compliance to the OSHA Bloodborne Pathogen Standard
- **Containment**—Seams comply with OSHA 3130 pertaining to containment of fluids
- **Leak Proof, Heat Sealed Seams**—Prevent leakage and minimize exposure to bloodborne pathogens
- **Resistant to Blood, Fluids, Body Fats and Chemicals**—For technician safety

### Eco-Friendly—Minimizing Chemicals of Concern and Waste

- Environment friendly olefin material construction is chlorine free (PVC-free)
- Low VOC emissions during incineration, for crematorium compatibility