

# EnviroTough-Bag™

## Eco-Friendly Heavy Duty Disaster Bags

**LDI**medical™  
Sustainable Technology · Clinical Innovation · Ultimate Performance

PART #	SIZE	PACKAGING
ETH2-3894	38" x 94" curved zipper with handles	6 bags per case



### Product Details—Our Strongest Body Bags

- **Superior Tensile Strength**—20 mil thickness, our strongest body bag for use in extreme situations
- **Heavy Duty, Rust Proof Zipper**—With two extra long pulls for easy gripping
- **Handles**—Six, integral leak proof handles with polypropylene webbing
- **Envelope Style Opening**—For easier insertion and removal of cadaver
- **20% Lighter**—Yet stronger than vinyl body bags—less weight, better performance

### Manage Risk and Comply with OSHA Standards

- **Compliance**—Aides compliance to the OSHA Bloodborne Pathogen Standard
- **Containment**—Seams comply with OSHA 3130 pertaining to containment of fluids
- **Meets CDC Minimum Thickness**—Guidelines for containment of Ebola
- **Leak Proof, Heat Sealed Seams**—Prevent leakage and minimize exposure to bloodborne pathogens
- **Resistant to Blood, Fluids, Body Fats and Chemicals**—For technician safety

### Eco-Friendly—Minimizing Chemicals of Concern and Waste

- Low VOC emissions during incineration, for crematorium compatibility
- Environment friendly olefin material construction is non-chlorinated (PVC-free)