

EnviroSorb[®] UOC

Reusable Underpads with Urine Odor Control

PART #	SIZE	PACKAGING
UOC-3436	34" x 36"	24 pads per case



*Urine Odor Control

Product Details–Discrete Control of Urine Odor

- **Rapidly Minimizes Urine Odor**–Proprietary material absorbs odor for patient and staff comfort
- **Reusable**–Launder up to 30 times with either bleach or peroxides
- **Cutting Edge Moisture-Barrier Technology**–Traps fluid for patient comfort and for leakage prevention
- **Fluid Absorbent Core**–Draws moisture away from patient
- **Soft, Fast Wicking Face Material**–Allows perspiration and moisture to escape
- **Lay Flat Technology**–No bunching or rolling for patient comfort
- **Unique Urine Neutralizing Chemistry** - Continues to work even after multiple launderings with bleach or peroxides.

Eco-Friendly Construction

- **Compliant**–Neutralizing chemical agent is Prop 65 compliant
- **Low Profile**–Bonded construction thins the pad with maximum absorption
- **Quilted**–To a blended poly/rayon soaker with a green film fluid barrier
- **Heavyweight**–55/45 cotton/polyester Twill